



950 N. Cole Rd. • Boise, ID 83704
4301 Garrity Blvd Suite 203 Nampa, ID 83687

Phone: (208) 323-CPAP • Fax: (208) 323-1615 • www.everythingcpap.com

How to CLEAN your PAP supplies

Daily: Clean your cushion with mask wipes, or warm soapy water and allow to air dry.

Weekly: Clean your mask frame, tubing, and water chamber with warm soapy water, allow to air dry.

Monthly: Sanitize your machine tubing and water chamber using admixture 3:1 ratio, 3 parts water 1 part white vinegar. Soak for an hour, thoroughly rinse, and allow to air dry.

How often to REPLACE your supplies

Every two weeks: You will need to replace your machine's disposable filter and nasal mask cushion, if applicable, every two weeks.

Monthly: You will need to replace your full face cushion monthly.

Every 3 months: You will need to replace your tubing and mask frame every 3 months.

Every 6 months: You will need to replace your water chamber and mask headgear every 6 months.

EverythingCPAP will call every 90 days to offer you renewable supplies. If you miss our call, please call us back to place your supply order.